


Konformitätserklärung

Certificate of Compliance

Körbchen nach USP/EP

Basket stirring element according to USP/EP

	USP/EP	ERWEKA
A	20,2±1,0 mm	20,2±1,0 mm
B	5,1±0,5 mm	5,1±0,5 mm
C	2,0±0,5 mm	2,0±0,5 mm
D	9,4 - 10,1 mm	9,4 - 10,1
E	25,0±3,0 mm	25,0± 1,0 mm
F	20,2±1,0 mm	20,2± 1,0 mm
G	22,2±1,0 mm	22,2±1,0 mm
J	37,0±3,0 mm	36,8± 1,0 mm
K	27,0±1,0 mm	27,0± 1,0 mm
L	400±40 µm	381±20 µm
M	0,22 - 0,31 mm	0,22 - 0,31 mm



Werkstoffe:	USP/EP	ERWEKA
Material:		
Spannfeder:	type 316	St.1.4571 ^①
Retention spring:		
Gewebe:	type 316	St.1.4401 ^②
Screen:		
Alle anderen Teile:	type 316	St.1.4571 ^①
All others:		

Für das Körbchen mit Schaft bescheinigen wir in alleiniger Verantwortung die Übereinstimmung mit den Forderungen der USP42, General Chapter <711> und der European Pharmacopoeia 10th Edition, Kapitel 2.9.3.

For the basket stirring element with shaft we declare in our sole responsibility the conformity with USP42, General Chapter <711> and the European Pharmacopoeia 10th Edition, Chapter 2.9.3.

Langen, 2020-05-27

i.A.


Günter Annacker
- Quality department -

^① German stainless steel, similar to AISI 316Ti

^② German stainless steel, similar to AISI 316